# Dentflex

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# Summary

|  |         |
|--|---------|
| 1. INTRODUCTION.....   | 33      |
| 2. DESCRIPTION AND PURPOSES.....   | 33      |
| 3. INSTALLATION.....   | 34      |
| 4. TECHNICAL FEATURES.....   | 35      |
| 4.1 Manufacturer's Declaration - Electromagnetic Emissions.....          | 36      |
| 4.2 Manufacturer's Declaration - Electromagnetic Emissions.....          | 37      |
| 4.3 Manufacturer's Declaration - Electromagnetic Immunity .....          | 38      |
| 4.4 Recommended Separation Distance .....                                | 39      |
| 4.5 Used symbology .....   | 40      |
| 4.5.1 Symbols in this Handbook Instruction.....                          | 40      |
| 4.5.2 Symbols on Packaging, on Appliance, and on Parts of Appliance..... | 40      |
| 5. HANDBOOK INSTRUCTIONS.....  | 41      |
| 5.1 Content.....   | 41      |
| 5.2 Preparation.....   | 42      |
| 5.3 Equipment Installation .....   | 42      |
| 5.3.1 Control Unit .....   | 42      |
| 5.3.2 Power supply.....  | 43      |
| 5.3.3 Command Pedal .....  | 43      |
| 5.3.4 Electric Micromotor .....  | 45      |
| 5.3.5 Holder .....   | 46      |
| 5.3.6 Peristaltic Pump.....  | 46      |
| 5.3.7 Irrigation Kit.....  | 47      |
| 5.3.8 Serum Stem .....   | 49      |
| 6. EQUIPMENT OPERATION.....  | 50      |
| 6.1 Display.....   | 50      |
| 6.1.1 Program .....  | 51      |
| 6.1.2 Gear .....   | 51      |
| 6.1.3 Speed.....   | 51      |
| 6.1.4 Torque .....   | 51      |
| 6.1.5 Pump.....  | 52      |
| 6.1.6 Rev .....  | 52      |
| 7. ENVIRONMENTAL PROTECTION .....  | 53      |
| 8. BIOCOMPATIBILITY .....  | 53      |
| 9. STORAGE .....   | 53      |
| 10. WARNINGS / CAUTIONS .....  | 53      |
| 11. MAINTENANCE / CLEANING / DISINFECTION / STERILIZATION.....           | 55      |
| 11.1 Maintenance.....  | 55      |
| 11.2 Cleaning, Disinfection and Sterilization .....                      | 55      |
| 11.2.1 Cleaning and Disinfection .....                                   | 56      |
| 11.2.2 Sterilization .....   | 56      |
| 12. FUSES' REPLACEMENT .....   | 56      |
| 13. POSSÍVEIS DEFEITOS / CAUSAS / SOLUÇÕES.....                          | 57      |
| 14. WARRANTY .....   | 57      |
| 15. DENTFLEX TECHNICAL SUPPORT CENTER.....                               | 58      |
| INSTRUCCIONES DE USO .....   | 59 - 86 |

## 1. INTRODUCTION

Congratulations on your purchase of the DForce 1000.

If you have any questions or problems while reading this handbook instructions, please contact Dentflex. Keep this manual in a safe place.

Dentflex reserves the right to incorporate, at any time and without prior notice, changes in the information and data indicated in this document. This handbook instructions has been developed with absolute rigor. However, even with our efforts, it is impossible to eliminate the possibility of inaccuracies. We are always grateful for any recommendations for improvement, so please contact Dentflex directly.

## 2. DESCRIPTION AND PURPOSES

This appliance is intended exclusively for dental applications. The DForce 1000 motor is designed for dentists' use only, together with medical instruments, it is used in implant surgery (perforation and implant insertion), periodontal surgery, dental wear and polishing. This equipment has a fundamental role in the patient's oral rehabilitation and oral health.

The DForce 1000 is a peripheral equipment that uses ISO 3964 coupling type surgical or dental handpieces. It has an electric micromotor, a peristaltic pump and a multifunction pedal to perform surgical and dental procedures.

### WARNING

- DForce 1000 is only intended for dental treatment in dental medicine. Its use for any other purposes is not permitted. It is dangerous to use it inappropriately.
- DForce 1000 should only be used by qualified dental professionals.

The equipment comes standard with 10 pre-configured procedures, such as: handpiece to be used, electric micromotor speed, peristaltic pump enabled or disabled and pump flow (amount of irrigated liquid). All functions are configured by means of a keypad and a display on the front of the equipment. The activation / deactivation of the electric micromotor and the peristaltic pump is done by foot control.

## WARNING

- Before operating the equipment, the user must configure it with its basic functions such as: handpiece to be used (gear), micromotor speed (speed) and working torque (torque).
- The equipment can be configured for up to 10 different procedures (program); after its configuration, it is operated by the pedal only.

The equipment functions are:

- Bone perforation performed by the electric micromotor in conjunction with a dental handpiece and drill;
- Operative field irrigation performed by the peristaltic pump;
- Implant's insertion by the electric micromotor in conjunction with a contra angle and an implant toolkit.

## NOTE

The DForce 1000 equipment is not supplied with handpieces, not even with an implant toolkit. There are no contraindications for the DForce 1000 equipment's uses.

## 3. INSTALLATION

- Storage and installation must be done in a place free from bad weather conditions.
- Make sure that the power grid where the engine will be installed is in accordance with the equipment's power supply according to the Technical Characteristics.
- Install the DForce 1000 on a flat surface, at an ergonomically suitable height for the operator. Allow enough space for handling accessories.
- The DForce 1000 supply voltage is 127 / 220V, and the voltage selection is manual, using selector switch.
- Position the cables so that they are free and out of places where they can be trampled and do not place any furniture on them, do not leave cables twisted or bent.
- Before connecting the power cable to an outlet, make sure that the DForce 1000 "on/off" switch is set to 0 (Off). Always connect the power cable to a to a three-pin bipolar outlet "2P + T".
- For this product's installation, please observe the Brazilian Technical Standard NBR 13.534 requirements: Electrical Installations for Health Care Establishments Safety Requirements.

## 4. TECHNICAL FEATURES

| Equipment Features   |  |
|--|--|
| Supply Voltage   | 127 / 220 V  |
| Input Power  | 120 VA   |
| Supply Voltage Selection   | Voltage selector switch  |
| Power Supply Frequency   | 60 Hz  |
| Dimensions   | 21,5 x 22,2 x 13,5   |
| Weight (approx. without accessories)                                       | 3,9 Kg   |
| Fuses  | 3A x 250V ~ / 5x20mm<br>(fast action - Type F)<br>breaking capacity<br>35A / 10 x In<br>(whichever is greater) |
| Equipment Rating   |  |
| Electric Shock Protection Type   | Class I  |
| Applied Part Protection Degree   | Type B   |
| Used in the presence of a flammable mixture such as air, oxygen or nitrous | equipment is not intended for use in these environments  |
| Equipment Operation Mode   | Non Continuous.<br>1 Minute (On) / 9 Minutes (Off)   |
| Mobility   | Portable   |
| Protection Degree Against Harmful Water Penetration                        | IP 31  |
| Environmental Conditions   |  |
| Storage  |  |
| Ambient Temperature  | from 5°C to + 50°C   |
| Relative Humidity  | from 10% to 80%  |
| Atmospheric Pressure   | from 500 to 1060 kPa<br>(from 375 to 795 mmHg)   |
| Operation  |  |
| Ambient Temperature  | from 10°C to 40°C  |
| Relative Humidity  | from 30% to 75%  |
| Atmospheric Pressure   | from 700 to 1060 kPa<br>(from 525 to 795 mmHg)   |
| Transportation   |  |
| Maximum Stacking   | 4 boxes  |

### NOTE

In case of interpretation doubts about any available data, contact our technical assistance by mail: [assistencia@dentflex.com.br](mailto:assistencia@dentflex.com.br).

## 4.1 Manufacturer's Declaration - Electromagnetic Emissions

**The Dforce 1000 is intended for use in the electromagnetic environment specified below.  
The DForce 1000 user must ensure that it is used in such environment.**

| Immunity Assays   | ABNT NBR IEC 60601<br>Test level                 | Compliance Level                                 | Electromagnetic<br>Environment – Guidelines   |
|---|--|--|---|
| Electrostatic discharge<br>(ESD) IEC 61000-4-2  | ± 2 KV per contact                               | ± 2 KV per contact                               | Floors should be made of wood, concrete or ceramic. If floors are covered with synthetic material, the relative humidity should be at least 30%   |
|   | ± 4 KV per contact                               | ± 4 KV per contact                               |   |
|   | ± 6 KV per contact                               | ± 6 KV per contact                               |   |
|   | ± 2 KV through the air                           | ± 2 KV through the air                           |   |
|   | ± 4 KV through the air                           | ± 4 KV through the air                           |   |
|   | ± 8 KV through the air                           | ± 8 KV through the air                           |   |
| Fast electric transients /<br>Burst Pulse Train<br>IEC 61000-4-4  | ± 2 KV on power lines                            | ± 2 KV on power lines                            | Energy supply quality should be that on from a hospital or commercial environment. If the DForce 1000 user requires continued operation during power interruption, is recommended that the DForce 1000 be powered by an uninterruptible power supply or by a battery. |
|   | ± 1 KV on input / output lines                   | ± 1 KV on input / output lines                   |   |
| Outbreaks IEC 61000-4-5   | ± 0,5 KV phase-phase                             | ± 0,5 KV phase-phase                             |   |
|   | ± 1 KV phase-phase                               | ± 1 KV phase-phase                               |   |
|   | ± 0,5 KV phase1 - GND                            | ± 0,5 KV phase1 - GND                            |   |
|   | ± 1 KV phase1 - GND                              | ± 1 KV phase1 - GND                              |   |
|   | ± 2 KV phase2 - GND                              | ± 2 KV phase2 - GND                              |   |
|   | ± 0,5 KV phase2 - GND                            | ± 0,5 KV phase2 - GND                            |   |
|   | ± 1 KV phase1 - GND                              | ± 1 KV phase1 - GND                              |   |
| Voltage drops, short<br>interruptions and voltage<br>variations on power supply<br>lines IEC 61000-4-11 | <5% Ut (> 95% voltage drop in Ut) for 0.5 cycle. | <5% Ut (> 95% voltage drop in Ut) for 0.5 cycle. | Energy supply quality should be that on from a hospital or commercial environment. If the DForce 1000 user requires continued operation during power interruption, is recommended that the DForce 1000 be powered by an uninterruptible power supply or by a battery. |
|   | 40% Ut (60% voltage drop in Ut) for 5 cycles     | 40% Ut (60% voltage drop in Ut) for 5 cycles     |   |
|   | 70% Ut (30% voltage drop in Ut) for 25 cycles..  | 70% Ut (30% voltage drop in Ut) for 25 cycles    |   |
|   | <5% Ut (> 95% voltage drop in Ut) for 5 seconds. | <5% Ut (> 95% voltage drop in Ut) for 5 seconds. |   |
| Magnetic field at power frequency (50/60 Hz)<br>IEC 61000-4-8   | 3 A/m  | 3 A/m  | Magnetic fields in feeding's frequency should be in typical hospital or commercial environment characteristic levels. Los campos magnéticos .   |

NOTA - Ut is the Alternating Current supply voltage before the test level application

## 4.2 Manufacturer's Declaration - Electromagnetic Emissions

The DForce 1000 is intended for use in the electromagnetic environment specified below.  
The DForce 1000 user must ensure that it is used in such environment.

| Emission Tests   | Conformity | Electromagnetic Environment – Guidelines  |
|--|------------|---|
| RF Emissions<br>ABNT NBR IEC CISPR11                                 | Group 1    | The Dforce 1000 FIX uses radio frequency (RF) energy only for its internal functions. However, their RF emissions are very low and are not likely to cause any interference to nearby electronic equipment. |
| RF Emissions<br>ABNT NBR IEC CISPR11                                 | Class B    |   |
| Harmonic Emissions<br>IEC 61000-3-2                                  | Class A    | DForce 1000 is suitable for use in non-residential establishments and those directly connected to the low voltage electricity public distribution network, which feeds buildings for domestic purposes.     |
| Emissions due to voltage fluctuation / scintillation<br>IEC61000-3-3 | According  |   |

## 4.3 Manufacturer's Declaration - Electromagnetic Immunity

DForce 1000 is intended for use in the electromagnetic environment specified below.  
The DForce 1000 user must ensure that it will be used in such environment.

| Immunity Assays   | ABNT NBR IEC 60601 Test level   | Compliance Level | Electromagnetic Environment – Guidelines  |
|---|---------------------------------|------------------|---|
| RF Conducted<br>IEC 61000-4-6   | 3 Vrms<br>0.15 KHz up to 80 MHz | [3] V            | Portable and mobile RF communication equipment should not be used near any part of the DForce 1000, including cables, with separation distance less than the recommended one, calculated from the equation applicable to the transmitter frequency. |
| <b>Recommended Separation Distance</b>  |                                 |                  |   |
| $d = [3,5/V1] \cdot \sqrt{P}$ $d = [3,5/3] \cdot \sqrt{P} = 1,16 \cdot \sqrt{P}$  |                                 |                  |   |
| $d = [3,5/E1] \cdot \sqrt{P}$ 80 Mhz up to 800 Mhz $d = [3,5/3] \cdot \sqrt{P}$ 80 Mhz up to 800 Mhz = 1,16 $\cdot \sqrt{P}$  |                                 |                  |   |
| $d = [7/E1] \cdot \sqrt{P}$ 800 Mhz up to 2,5ghz $d = [7/3] \cdot \sqrt{P}$ 800 Mhz up to 2,5ghz = 2,33 $\cdot \sqrt{P}$  |                                 |                  |   |
| where P is the transmitter's maximum rated output power in watts (W), according to the transmitter's manufacturer, and d is the recommended separation distance in meters (m).  |                                 |                  |   |
| It is recommended that the field strength established by the RF transmitter, as determined by an on-site electromagnetic inspection, should be less than the compliance level in each frequency range <sup>b</sup> .  |                                 |                  |   |
| Interference may occur around the equipment marked with the following symbol:   |                                 |                  |   |
|   |                                 |                  |   |
| NOTE 1 - For 80 MHz to 200 MHz - 200 MHz to 1000 MHz - 1 GHz to 2.5 GHz, the higher frequency range is applied.   |                                 |                  |   |
| NOTE 2 - These guidelines may not apply in all situations. Electromagnetic propagation is affected by the absorption and reflection of structures, objects and people.  |                                 |                  |   |
| <ul style="list-style-type: none"> <li>• The field strengths established by the fixed transmitters, such as: base stations, telephone (cellular / wireless) and terrestrial mobile radios, amateur radio, AM and FM radio transmission and TV transmission; cannot be predicted theoretically with precision. To evaluate the electromagnetic environment due to fixed RF transmitters, an electromagnetic site inspection is recommended. If the field strength measurement at the location where the DForce 1000 is used exceeds the compliance level used above, the DForce 1000 should be observed to verify if the operation is normal. If abnormal performance is observed, additional procedures may be required; such as reorienting or replacing the DForce 1000.</li> </ul> |                                 |                  |   |
| <ul style="list-style-type: none"> <li>• Frequency range above of 150 kHz to 80 Mhz, the field strength should be less than [3] V / m.</li> </ul>   |                                 |                  |   |

## 4.4 Recommended Separation Distance between portable and mobile RF communication equipment and the DForce 1000

The DForce1000 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The implant motor user can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communication equipment (transmitters) and the DForce 1000 as recommended below, according to the communication equipment's maximum output power.

| Transmitter's maximum rated output power<br>W | Separation distance according to transmitter frequency<br>m         |   |   |
|---|---|---|---|
|   | 150 kHz up to 80 MHz<br>$d = \left[ \frac{3,5}{3} \right] \sqrt{P}$ | 80 Mhz up to 800 Mhz<br>$d = \left[ \frac{3,5}{3} \right] \sqrt{P}$ | 80 Mhz up to 2,5 Ghz<br>$d = \left[ \frac{7}{3} \right] \sqrt{P}$ |
| 0,01  | 0,11  | 0,11  | 0,23  |
| 0,1   | 0,36  | 0,36  | 0,73  |
| 1   | 1,16  | 1,16  | 2,23  |
| 10  | 3,68  | 3,68  | 7,37  |
| 100   | 11,66   | 11,66   | 23,33   |

For transmitters with a maximum rated output power not listed above, the recommended separation distance d in meters (m) can be determined by the equation applicable to the transmitter frequency, where P is the transmitter's maximum rated output power in watts (W) according to the transmitter's manufacturer.

NOTE 1 - At 80 Mhz and 800 Mhz, the separation distance is applied to the highest frequency range.


NOTE 2 - These guidelines may not apply in all situations. Electromagnetic propagation is affected by the absorption and reflection of structures, objects and people.

## 4.5 Simbologia Utilizada

### 4.5.1 Symbols in this Handbook Instruction

|   |  |  |  |
|---|--|--|--|
| <br><b>WARNING</b> | If the instructions given in the WARNINGS are not followed correctly, this may cause equipment failure and endanger user and patients. | <br><b>NOTE</b> | NOTES contain supplementary information and explanations about the equipment's use and operation |
|---|--|--|--|

### 4.5.2 Symbols on Packaging, on Appliance, and on Parts of Appliance

|  |  |   |  |
|--|--|---|--|
|    | General Warning Symbol   |    | General Mandatory Action Symbol  |
|    | General Prohibition Symbol   |    | See Handbook Instruction   |
|    | Manufacture Date:<br>Indicates the equipment's manufacture date  |    | Fuse   |
|    | Equipment On   |    | Equipment Off  |
|    | Type B applied part: Indicates the protection degree against electric shock from the part of the equipment that comes in contact with the patient            |   |  |
|   | Autoclave Sterilization Condition  |    | Maximum Sterilization Temperature (attention to the temperature indicated by the symbol) |
|   | Temperature Limits: Indicates the temperature limits for packaging's transportation and storage  |   | Fragile: The contents in this package are fragile and must be transported with care      |
|    | Keep Out of the Rain: Packing must not be carried in the rain  |  | This Side Up: Indicates the correct position for packaging's transportation              |
|    | Maximum Stacking: Maximum number of identical packages that can be stacked one on top of the other. For this equipment, the stacking limit number is 4 units |   |  |
|    | Indicates that the product should be taken to a special waste disposal site at the end of its useful life. Applies to: device and accessories                |   |  |

## 5. HANDBOOK INSTRUCTIONS

When opening the package and before installation, check the product for damage or missing parts. If shipping damage or missing parts are detected, a claim must be made to the distributor within 24 hours.


### WARNING

- Do not place the appliance in humid places or where it may be in constant contact with liquids.

### 5.1 Content

The DForce 1000 is composed by:

## NOTE

The handbook instructions and the warranty statement are included with the product

## 5.2 Preparation

1. Carefully remove the unit and other parts from the package and place everything on a flat surface.
2. Verify that all components described in CONTENT are present.
3. Verify that the serial numbers on the product, package, and warranty statement are identical to each other.
4. Environmental conditions for operation: indoor application, ambient temperature from 15°C to 40°C and relative air humidity from 30% to 75%.
5. If stored under ambient conditions of 5°C to 50°C and relative humidity from 10% to 80%, the original package materials may be reused in the future.

### WARNING

- Enclosure protection rating is IP 31.

## 5.3 Equipment Installation

### 5.3.1 Control Unit

The Control Unit was designed to control the electric micromotor in surgical and dental procedures, allowing torque control, speed control, reduction control, micromotor reversal and pump flow control.

For the production of this equipment the following materials were used: Stainless Steel, Aluminum, Polyurethane Box and Electro - Electronic Components. (Fig. 01)




Fig. 01

### 5.3.2 Power supply

Motor operation is by AC power only. (Fig. 02)

a. Plug the power cord into a power outlet and insert the connector into the back of the machine.

b. Press the button on the back to turn the unit on or off.

PP Power Cord, 3-Way x 1mm, Length 2.3 meters. It has 2 plugs of 3 contacts being one (M) male and the other (F) female. (Fig. 03)

This cable makes electrical connections and disconnections (PHASE, NEUTRAL and GROUND) between the equipment and the electrical installation point where it will be connected.

#### WARNING

- Before plugging in the power cord, make sure that the machine's power switch on / off is set to 0.
- To disconnect, always pull on the center of the connector. Never pull on the cable.
- If you have a problem, contact the Technical Support Center.

### 5.3.3 Command Pedal

Ergonomic pedal consisting basically by 1 foot pedal and 2 buttons, in blue and yellow colors, fixed in an aluminum and carbon steel base. The command pedal is connected to the central control box via a connector.

#### WARNING

- The command pedal has IPX1 protection degree and cable of 2,30 meters in length.
- In the event of a power outage or breakdown, it is necessary to release the foot pedal and turn it on again in order to return the equipment to normal operation.
- To ensure a long service life, avoid bending the cables.
- Do not continually push the pedal cord against firm objects.

For the connection between the pedal and the control unit, attach the pedal connector (a) to the motor connection (b) located at the rear of the motor. (Fig. 04)


Fig. 04

### c. Speed Control

The pedal has the analogically function of controlling the micromotor rotation; pressing it will gradually accelerate according to the range of rotation, which can be controlled by the professional according to the pressure exerted on the pedal. The applied rotation will be shown on the display and its rotation range will vary according to the used reduction.

### d. Yellow Button

#### Peristaltic Pump Activation / Deactivation

By holding this button down for approximately 4 seconds, the peristaltic pump will be enabled and at the display, below the FLOW function, will have an indicator signaling its activation. By repeating this procedure, the pump will be disabled, and the indicator at display will disappear.

#### Program Changes

DForce 1000 has 10 programs that can be parameterized by the user as needed (record torque, speed, flow and reduction). Whenever the professional is performing an operation that has already been parameterized in one of the programs, simply search for the desired option by touching the yellow button

### e. Blue Button

#### Anticlockwise Activation

Holding down this button will trigger the electric micromotor counterclockwise, and you can see its direction of rotation on the display. Prior to starting the micromotor, the equipment beeps 5 times to alert the operator that reversal will occur. When the foot is removed from the button, the motor automatically stops turning and the rotation's direction returns to the clockwise movement and is also indicated on the display.

### 5.3.4 Electric Micromotor

Insert the micromotor cable connector into the 8-pin metal socket on the front of the unit with the arrows facing up. (Fig. 05)


Fig. 05

The Electric Micromotor operates at a maximum speed of 35000 RPM, 30V working voltage, has an ISO 3964 coupling and a 2 meter long cable. The user must select, through the GEAR button, the reduction to be used in the handpiece coupled to the electric micromotor, which can be: 1: 1; 16: 1; 20: 1; 32: 1; 1: 2.

#### WARNING


- If the professional selects a reduction and fits a handpiece with a different reduction than the one selected, the equipment will not perform as expected and may be damaged.

#### NOTE

The handpiece to be used must fit to the front face of the Electric Micromotor, thus preventing premature wear of the shaft coupling.

Speed increase or reduction depends on the attached handpiece.

- The Electric Micromotor must always be sterilized in autoclave after being used at 135 °C, pressure of 2.2 bar and cycle of 35 minutes.

- Never lubricate or allow lubricating oil to drain into the Electric Micromotor, such a procedure may cause irreversible damage to the instrument.

#### WARNING

- The Electric Micromotor has Type B applied part, because it does not come into direct contact with the patient; the products that come in contact with the patient are: the handpiece and the drill.
- The Electric Micromotor protection degree is IPX3.

### 5.3.5 Holder


Fig. 06

Place the handpiece correctly connected at the resting Electric Micromotor in the holder. (Fig. 06)

#### NOTE

The holder serves as rest support to the Electric Micromotor and to the handpiece, thus preventing surface scratches and falls.

This holder should not be autoclaved.

### 5.3.6 Peristaltic Pump

The Peristaltic Pump performs the saline's pumping, whose function is to perform the refrigeration of the place to be worked. It is composed by metal alloys, such as: anodized aluminum and chrome brass. Its drive is done through an electric motor coupled inside the central control box.


Fig. 07


#### WARNING

- The peristaltic pump should only be handled when it is stopped, if necessary, first turn off the equipment.

### 5.3.7 Irrigation Kit

Indicated for use in the Dentflex peristaltic pump for irrigation in implant procedures and surgeries. Is composed by:

**a. 0.5 m silicone hose**

This tube goes through the peristaltic pump, at its ends are connected: the serum bag connector and the irrigation reducer.

**b. Serum Bag Connector**

Its function is to connect the peristaltic pump inlet to the saline reservoir.

Observe the arrows on the bottom of the cassette before connecting the serum bag connector.

**NOTE**

Before piercing the serum bag with the connector, make sure the irrigation kit is assembled correctly as shown above. With the kit assembled, pierce the serum bag with the perforation points facing up, then hang it on the serum stem.

**c. Irrigation reducer**

Its function is to connect the peristaltic pump output (0.5 m tube) to the 2.5 m tube.

**d. Containment nipple**

Its function is to prevent the hose from moving during the peristaltic pump operation.

**e. 2.5 m silicone hose**

This tube is connected to the irrigation reducer and then to the handpiece, bringing saline from the peristaltic pump outlet to the handpiece that will be used.


**WARNING**

- If the Contra Angle has internal irrigation, the tube is connected to the irrigation reducer and subsequently to the "Y" connector, bringing the saline from the peristaltic pump outlet to the handpiece that will be used.


## Assembly and Disassembly

- 1 - Lock Nut
- 2 - Peristaltic Pump Arm
- 3 - Pump Reel
- 4 - Fitting Clamp

### Irrigation Kit Assembly (Fig. 08):

- Turn Lock Nut (1) counterclockwise;
- Open the Peristaltic Pump Arm (2);
- Insert the 0.5 m silicone hose into the Fitting Clamp (4) with the Containment Nipple to the left;
- Pass the 0.5 m silicone hose over the Pump Reel (3);
- Lower the Peristaltic Pump Arm (2);
- Screw in the Lock Nut (1), turn it clockwise;
- Insert the 2.5 m silicone hose at the tip of the Irrigation Reducer.

### Irrigation Kit Disassembly (Fig. 09):

- Turn Lock Nut (1) counterclockwise;
- Open the Peristaltic Pump Arm (2);
- Remove the 0.5 m silicone hose from the Fitting Clamp (4) and from above the Pump Reel (3);
- Lower Peristaltic Pump Arm (2);
- Screw in the Lock Nut (1), turn it clockwise.


Fig. 09

### WARNING

- All Irrigation Kit components are autoclavable at 135°C, 2.2 bar pressure and 35 minute cycle.

### 5.3.8 Serum Stem

The Serum Stem should be used to fix the saline bag. It is made of chromed stainless steel. The support located at the rear of the unit serves to secure the stem. (Fig. 09) The Serum Stem should always face outwards. (Fig. 10). If a leak occurs in the saline bag, this positioning prevents fluid keep in contacting with the motor.


Fig. 09


Fig. 10

### WARNING

- If the motor keep in touch with any kind of liquid and the device has a problem, contact the Technical Support Center.
- Improper use will void the warranty.


## 6. EQUIPMENT OPERATION

After installing the equipment correctly press the Power key On / Off (A) located on the rear of the motor.


On / Off

Used to turn the appliance on or off

### 6.1 Display

#### KEYBOARD

1. Program - Programming selection key
2. Gear - Handpiece selection Key
3. Speed - Speed selection key
4. Torque - Torque selection Key
5. Pump - Peristaltic Pump adjustment Key
6. Rev - Rotation direction selection key

#### GRAPHIC PANEL

- a. Program- Selected Program
- b. Gear - Selected Handpiece
- c. Speed - Selected Speed
- d. Torque - Selected Torque
- e. Pump - Drive and Flow
- f. Engine - Selected Rotation

## 6.1.1 Program

DForce 1000 has 10 program options, which can be adjusted as required by the user.

To select the desired program click on the arrows.

The selected option can be viewed through the Graphic Panel.

Another way to change programs is through the yellow pedal button, as shown on page 15.

## 6.1.2 Gear

DForce 1000 has 05 options for handpieces' use that can be adjusted as required by the user:

- Handpiece Transmission: 1:1
- Reduction Handpiece: 16:1, 20:1 and 32:1
- Multiplier Handpiece: 1:2o:

To select the desired program click on the arrows.

The selected option can be viewed through the Graphic Panel.

- Only handpieces with these transmission, reduction and multiplication values can be connected to the DForce 1000.
- Never use handpieces with values other than the selected option. Before use check the mode that is indicated in the Graphic Panel.
- Improper use will void the warranty.

## 6.1.3 Speed

To select the desired program click on the arrows.

The selected option can be viewed through the Graphic Panel.

The micromotor operating speed value depends on the Gear option selection as shown in the following table:

| Gear | Speed (RPM)  | Variations Range (RPM)        |
|------|--------------|-------------------------------|
| 1:1  | 1000 – 35000 | 1000                          |
| 16:1 | 15 – 2000    | <100 RPM = 5<br>>100 RPM = 25 |
| 20:1 | 15 – 1750    |                               |
| 32:1 | 15 – 1000    |                               |
| 1:2  | 2000 - 70000 | 2000                          |

## 6.1.4 Torque

To select the desired program click on the arrows.

The selected option can be viewed through the Graphic Panel.

The torque value depends on the Gear option selection. When the Gear value of the handpiece attached to the micromotor is 16: 1, 20: 1 or 32: 1 reduction, the torque can be adjusted as shown in the following table:

| Gear | Speed (RPM)  | Torque (Motor) | Variations Range (N.cm) |
|------|--------------|----------------|-------------------------|
| 1:1  | 1000 - 35000 | 100%           | —                       |
| 16:1 | 15 - 2000    | 5 - 55 Ncm     | 5                       |
| 20:1 | 15 - 1750    | 5 - 55 Ncm     |                         |
| 32:1 | 15 - 1000    | 5 - 55 Ncm     |                         |
| 1:2  | 2000 - 70000 | 100%           | —                       |

### 6.1.5 Pump

To select the desired program click on the arrows.  
The selected option can be viewed through the Graphic Panel.

By pressing the FLOW key it is possible to adjust the pump flow rate with the following parameters:

| Minimum Flow | Max Flow     | Variations Range |
|--------------|--------------|------------------|
| 25 ml / min  | 150 ml / min | 25 ml / min      |

By pressing the ENABLE / DISABLE key it is possible to enable as well as disable the Pump.

#### NOTE

Another way to enable and disable the pump is by using the yellow pedal button as shown on page 15.

### 6.1.6 Rev

This key is used to change the electric micromotor rotation direction.

Pressing the REV key, the micromotor rotation will be counterclockwise and in the Graphic Panel you will see the rotation direction.

To return clockwise, press the REV key again. You will also be able to view the direction of rotation through the Graph Panel.

#### NOTE

Another way to activate the reversal rotation direction mode is by pressing the blue pedal button as shown on page 15.

#### WARNING

- When this button is pressed or the blue button on the pedal is pressed, the unit will beep to warn the operator that the Micromotor is counterclockwise.

## 7. ENVIRONMENTAL PROTECTION

Like any other electronic device, the DForce 1000 FIX is composed by several materials, such as: plastics, metals, rubbers, among others. For this reason, the product's definitive disposal should not be made in ordinary trash, but in an appropriate place. Find a competent body in your city to assist you.

## 8. BIOCOMPATIBILITY

The materials used in the handpiece are inert, do not cause allergic reactions.

## 9. STORAGE

Keep in a place protected from the sun, heat, dust and moisture; away from chemicals and cleaning agents. The equipment must be stored, transported in its respective carrying case, at ambient temperature and free of excessive pressure.

## 10. WARNINGS / CAUTIONS

### Warnings

1<sup>a</sup> Considering the number of radiofrequency sources increase, present in our daily life, the equipment's interference occurrence probability increases. In this way, it is always important to avoid that electronic products, especially medical equipment, are close to radiofrequency devices. It is, therefore, necessary that electronic devices are electromagnetically compatible, that is, they must be designed to withstand a certain limit of external interference and should not emit interference beyond the controlled limits. The electromagnetic compatibility standard NBR-IEC 60601-1-2 applied to electro-medical equipment prescribes the limits of the electromagnetic field strength that these devices must withstand and the limits that they can emit. This equipment complies with the requirements for electromagnetic compatibility In accordance with NBR-IEC 60601-1-2, CISPRII, Class B. Before installing the equipment, refer to the handbook instructions for 'Technical Characteristics' and 'Warnings / Cautions / Recommendations, and use the connection to an electrical network that is within the normative requirements and standards.

2<sup>a</sup> It is not advisable to place electro-medical equipment to work together with mobile or portable radiofrequency communication equipment.

3<sup>a</sup> The use of parts, accessories or spare parts that make up the product must be carried out according to Dentflex's instructions. The use of parts, accessories or spare parts not specified / approved by Dentflex may result in increased emissions or decreased immunity of the engine.

4<sup>a</sup> DForce 1000 must not be used adjacent or stacked to other equipment.

5<sup>a</sup> The use of accessories and cables other than those for which the equipment is designed may significantly degrade emissions and immunity performance.

6<sup>a</sup> All parts and accessories supplied with the equipment have been tested and approved by Dentflex. Using other parts, other than those described, may compromise the equipment's safety and performance.

7<sup>a</sup> The Motor DForce 1000 is not supplied with dental handpieces, not even drills and implant kits. In order to use these accessories, the user needs to be aware of their use and can acquire this knowledge in specific courses.

8<sup>a</sup> To avoid electric shock risk, this equipment must be connected to a grounded power supply only.

9<sup>a</sup> According to Standard 60601-1 p. 52 (item 7.9.3.1 1st par. trace 8).

10<sup>a</sup> No modifications to this equipment are allowed.

11<sup>a</sup> Do not modify this equipment without the manufacturer's explicit permission.

12<sup>a</sup> Even if authorized by the manufacturer, if this equipment is modified, appropriate inspections and tests must be carried out to ensure continued safety of use.

13<sup>a</sup> When using the motor, the operator should not touch the equipment's connectors.

14<sup>a</sup> In the event of a power outage or breakdown, it is necessary to release the foot pedal and turn it on again in order to return the equipment to normal operation.

## Cautions

- When positioning the equipment properly to work, leave enough space to unplug the power cord plug from the equipment's outlet;
- Power cord length is 2.3 meters;
- Do not allow any liquid to enter inside the equipment;
- Do not allow that the equipment fall to the ground;
- Do not allow oil to enter inside the Electric Micromotor, this implies loss of warranty;
- Avoid working with the overheated Electric Micromotor. If this happens, stop using it and use it alternately;
- Do not use any components not mentioned in this handbook instruction. Using different components may damage the equipment;
- Do not use any component of the equipment if it is damaged;
- The use of any part, accessory or unspecified material is prohibited;
- Never use worn drills;
- Never lubricate the Electric Micromotor;
- Never connect or remove a Contra Angle with the Electric Micromotor turned on;
- Never disconnect the ground wire removing it by the cord or using adapters. This procedure, if performed, may cause damage to the control unit which will not be covered by the warranty;
- When unplugging the power cord, use the plug. Never pull on the cord.

## 11. MAINTENANCE / CLEANING / DISINFECTION / STERILIZATION

DForce 1000 regular maintenance is a fundamental requirement for hygienic work with the device. Therefore, follow the maintenance, cleaning, disinfection and sterilization instructions in this chapter. As a manufacturer, Dentflex is responsible for the equipment's technical characteristics and safety only in cases where the unit has been used in accordance with the use instructions contained in this manual; where maintenance, repairs and modifications have been made by the factory or expressly authorized agents , and in which components that may cause safety and operating risks have been replaced by original spare parts in a fault event.

### 11.1 Maintenance

Inspect the micromotor and pedal charger cables at least once every 6 months. If any wear on the cable sheath is detected, contact the Dentflex Technical Support Center. Any lubricants or oils may contaminate the micromotor and cause unsafe operation. This will void the warranty. If you notice any liquid or smoke coming from the control unit, immediately disconnect the device from the battery charger and contact the Dentflex Technical Support Center.

#### NOTE

Repair and maintenance work should only be performed by the Dentflex Technical Support Center.

#### WARNING

- Prevent that the lubricant from the contra angle entering the micromotor. Eliminate excess lubricant by injecting compressed air into the back of the contra angle, blow approximately 5 seconds before connecting to the micromotor.
- The micromotor should never be lubricated.
- The warranty is void if any non-original parts are used. Example: micromotor, pedal etc.

### 11.2 Cleaning, Disinfection and Sterilization

As part of your responsibility, use only validated methods for cleaning, disinfection and sterilization. Ensure that disinfecting and sterilizing appliances are subject to regular technical inspection and that valid parameters are maintained during each cycle. Always follow the valid regulations and legal provisions that apply to your office or clinic.

#### NOTE

For your own safety when handling contaminated accessories always wear gloves, goggles and a mask.

## 11.2.1 Cleaning and Disinfection

The cables, the micromotor, the pedal and the surface of the appliance should be cleaned with a cleaning paper or a 70% alcohol dampened cloth. Soon after, use a damp cloth and finally a dry cloth. Perfect cleaning and disinfection are fundamental requirements for effective sterilization.

After cleaning / disinfection inspect all parts.

### WARNING

- Never immerse the equipment in ultrasonic cleaners.

## 11.2.2 Sterilization

The sterilization procedure applies only to the Electric Micromotor.

Follow sterilization procedures as per autoclave manufacturer.

The maximum sterilization temperature is 135°C.

### WARNING

- Never place the micromotor in ultrasonic cleaners. No components of the DForce 1000 should be sterilized (except the electric micromotor).
- Cables should never be autoclaved.
- The device control unit is not waterproof. Never apply liquids or sprays directly on the unit (especially on the dial) or near electrical connections.

## 12. FUSES' REPLACEMENT

The Fuse Holder is used to facilitate the replacement of a fuse that has been blown for any anomaly. DForce 1000 uses 5 / 20mm glass fuse, 250V voltage and fast acting 3rd current (Type F). The fuse holder is sent to the customer with two fuses, one of which is spare part. Low Breaking Capacity 35A / 10 x In (Rated Current) whichever is greater. Fuse's use not specified by manufacturer will void warranty.

|  |                                 |   |                                       |   |
|--|---------------------------------|---|---------------------------------------|---|
|  |                                 |   |                                       |   |
| Use a key to open the drawer until it is released. | Pull the drawer with the fuses. | Remove the damaged fuse from its housing. | Insert the new fuse into the housing. | Replace the drawer by pressing it until it locks. |

NOTE  
The fuse's characteristics can be found in the "Technical Features" section.

## 13. POSSIBLE DEFECTS / CAUSES / SOLUTIONS

If DForce 1000 does not work properly, first check the indications in the table below.  
If the problem persists, contact the Dentflex Technical Support Center.

| POSSIBLE DEFECTS   | CAUSES  | SOLUTIONS   |
|--|---|---|
| Equipment does not turn on   | Power cable improperly connected<br>Damaged power cable<br>Blown fuse<br>Power switch off       | Check connection<br>Replace damaged cable<br>Replace blown fuse<br>Turn on power switch   |
| Command pedal does not work  | Pedal cable improperly connected<br>Damaged component   | Check connection<br>Forward product to the Dentflex Technical Support Center  |
| Micromotor does not work   | Micromotor cable improperly connected<br>Damaged component                                      | Replace damaged cable<br>Forward product to the Dentflex Technical Support Center   |
| Contra angle does not fit  | Contra angle coupling is not in accordance with ISO 3964 standard                               | Use contra angles' coupling which is in accordance with ISO 3964 standard   |
| Peristaltic pump does not work   | Pump disabled<br>Pump locked due to the presence of air inside<br>Damaged component             | Turn on the pump through the pedal and pressing the yellow button or by the display through the "pump" button<br>Turn on the pump at maximum flow in order to eliminate the air<br>Forward product to the Dentflex Technical Support Center |
| Electric micromotor overheating  | Equipment's continuous use<br>Use contra angle with a different reduction than the selected one | Select the correct reduction  |
| Drill rotation does not match the rotation selected on the central control box | Use contra angle with a different reduction than the selected one                               | Select the correct reduction  |
| Insufficient torque  | Use contra angle with a different reduction than the selected one                               | Select the correct reduction  |

## 14. WARRANTY

The product manufactured by Dentflex is warranted, if it has been installed in accordance with the guidelines described in this manual. The warranty will be void and the manufacturer disclaims any liability in the following cases:

- If the user is not meeting the requirements set forth in the Use Instructions and has violated the equipment.
- Application of the appliance for any purpose other than its intended purpose as indicated in the instructions for use and maintenance.
- Connecting the motor to a power source not complying with IEC 364.
- Modifications and repairs not made by the Dentflex Technical Support Center.
- Use of non-original parts or components not described in chapter 9.1 Contents.
- Instrument fracture due to misapplication.
- Damage to accessories or appliance due to sterilization in violation of this manual.

For the warranty period's information, see the Product Certificate.

#### NOTE

Submit the Invoice to the Technical Assistance during the warranty period.

We assume no responsibility for repairs made outside of the Dentflex Technical Support Center.

#### WARNING

- Any maintenance or intervention on the equipment can only be carried out by the Dentflex Technical Support Center.

## 15. DENTFLEX TECHNICAL SUPPORT CENTER

With the purpose of offering a service that guarantees 100% satisfaction, we present to the customer the Dentflex Technical Support Center, an innovation in the service that came to provide agility, safety and transparency in the maintenance process of your Dentflex product.

We have an infrastructure and a qualified professionals' team, guaranteeing the most complete technical assistance for the entire Dentflex product line.

- Trained technicians and engineers;
- Last generation machinery and equipment;
- Original components.

To use the service offered by the Dentflex Technical Support Center, please, contact us:  
E-mail: [sac@dentflex.com.br](mailto:sac@dentflex.com.br) / Phones: + 55 (16) 9 9404-4435 and + 55 (16) 2133-0000.

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